



Product Code and Specification

Product Model	Specifications		
	Width (±5mm)	Length (±5mm)	Thickness (mm)
RDS-1	15	20	0.1-0.5
RDS-2	20	30	
RDS-3	30	40	
RDS-4	40	60	
RDS-5	60	60	
RDS-6	60	80	
RDS-7	60	140	
RDS-8	80	80	
RDS-9	80	120	
RDS-10	100	150	
RDS-11	150	150	


References

- Shi Z, Xu T, Yuan Y, Deng K, Liu M, Ke Y, et al. A new absorbable synthetic substitute with biomimetic design for dural tissue repair. *Artif Organs*. 2016;40:403-13.
- Novel Regenerative Nanofibrous Bio-device for Dural Defect Repair. Congress of Neurological Surgeons Annual Meeting, Washing DC, USA, 2011.
- In-vitro and clinical study on a novel synthetic absorbable biomimetic dural substitute. European Society for Pediatric Neurosurgery(ESPN) Congress, Rome, Italy, 2014.
- Electrospun Fibrous Mats with High Porosity as Potential Scaffolds for Skin Tissue Engineering. *Biomacromolecules*, 2008,9(7):1795-1801.
- Development of Novel Nanofibrous Dural Substitute for Dural Repair. The 14th World Federation of Neurological Societies Interim Meeting, Pernambuco, Brazil, 2011.
- Francesco Zenga , et al. Nanofibrous Synthetic Dural Patch for Skull Base Defects: Preliminary Experience for Reconstruction after Extended Endonasal Approaches. *Journal of Neurological Surgery Reports* 2016;77:e50-e55.
- Kunxue Deng, Xun Ye, Yaya Yang, Man Liu, et al. Evaluation of efficacy and biocompatibility of a new absorbable synthetic substitute as a dural onlay graft in a large animal model. *Neurological Research* 2016.


Medprin Biotech GmbH
 Address: Gutleutstraße 163-167, 60327 Frankfurt am Main, Germany
 Tel: +49-69-580059970 Fax: +49-69-580059971
 Email: Sales@medprin.com


ReDura™

Biomimetic-Synthetic-Absorbable
Dural Substitute



- Rapid Repair & Regeneration
- Long Term Safety
- Excellent Handling & Conformability
- High Strength & No-swelling
- Onlay & Suture






ReDura™


ReDura™ is made of poly-L-lactic acid(PLLA), which is FDA approved and has been applied worldwide in medical devices and proven its biocompatibility and non-toxicity.

Rapid repair and regeneration


ReDura™ resembles microstructure of human extracellular matrix(ECM) which provides ideal scaffold for cell proliferation and tissue regeneration.



Microstructure of native dura



Microstructure of ReDura™



Cells adhere tightly to the nanofibers of ReDura™

Prevent cerebrospinal fluid (CSF) leakage

ReDura™ is hydrophobic and acts as a physical barrier to prevent CSF leakage.

No liquid leakage using ReDura™



Excellent conformity

After hydration, ReDura™ becomes soft and conforms well to brain contour.



Native dura



ReDura™

High Strength


ReDura™ is of high mechanical strength.

Product	Strength (N/cm ²)
ReDura™	1.5-2
Control group	~0.45

Anti-adhesion


ReDura™: smooth surface of brain tissue with no adhesion to implanted material.

Control group: adhesion of native dura to brain tissue.



Clinical Guide


Suture procedure:



- 1 Trim ReDura™ into a suitable shape as needed.
- 2 Place ReDura™ onto the dural defect.
- 3 Fix ReDura™ into place by routine suturing with 4-0 sutures, or sutures of the surgeon's choice. During suturing, pinholes should stay 2-3mm away from the edge of ReDura™ to ensure a watertight closure.
- 4 Suturing Finished.

◆ ReDura can be sutured well without hydration, but it is highly recommended to hydrate it for optimal handling properties.

Onlay procedure:



- 1 Place the ReDura™ into normal or cold saline.
- 2 Press onto ReDura™ to ensure it absorbs the saline, remove it once it becomes translucent.
- 3 Trim ReDura™ into a suitable shape, the edges of ReDura™ should be beyond the defect area by 20-30mm.
- 4 Apply ReDura™ to cover the dural defect, ensuring the overlap of 20-30mm- a Valsalva manoeuvre can be performed to ensure a watertight closure. Sealant or sutures can be applied if a CSF leak is detected.