



POLYMER CLIP CLOSURE SYSTEM

Polymer Clip Closure System is designed to be implanted around blood vessels or a tissue bundle for permanent occlusion.

The ligation closure system is made of high-quality medical grade material to ensure patient safety.

A wide range of clip applicators are designed with colour coded systems for easy identification.

High-quality medical grade stainless steel is used on the reusable open surgery and endoscopic surgery applicators to ensure durability in the operating room and in the CSSD.



POLYMER CLIP CLOSURE SYSTEM AND APPLICATORS

Polymer ligating clips are made of non-absorbable Polymer, inert and are non-toxic. Clips have a high level of biocompatibility along with the structural and dimensional stability over time, necessary for long-lasting patient safety.

Designed with a self-locking catch mechanism and Integrated teeth to ensure stability and prevent slippage on the tissue. Tactile feedback confirms jaw seating and secure vessel placement. Precision bosses give security in the applicator which prevent dropouts. Non-conductive, radiolucent, and does not interfere with CT, MRI, or X-ray diagnostics.

Available clip sizes: Medium-Large, Large and Extra-Large. Clips are contained in colour coded cartridges for ease of identification.

Applicators for open and endoscopic surgery equipped with precision anchoring system that allows clips to remain stable without risk of accidental loss.

ORDERING CODES

DESCRIPTION

Polymer Clips

IM-KM-POLML
IM-KM-POLLG
IM-KM-POLXL

Polymer Ligation Clip Medium/Large Green
Polymer Ligation Clip Large Purple
Polymer Ligation Clip Extra Large Gold

Polymer Clip Applicators Laparoscopic Surgery

IM-KM-CA.ML.P33
IM-KM-CA.L.P33
IM-KM-CA.XL.P33
IM-KM-CA.ML.P44
IM-KM-CA.L.P44
IM-KM-CA.XL.P44

Polymer Laparoscopic 33cm Clip Applicator Medium/Large Green
Polymer Laparoscopic 33cm Clip Applicator Large Purple
Polymer Laparoscopic 33cm Clip Applicator Extra Large Gold
Polymer Laparoscopic 44cm Clip Applicator Medium/Large Green
Polymer Laparoscopic 44cm Clip Applicator Large Purple
Polymer Laparoscopic 44cm Clip Applicator Extra Large Gold