

WHEN YOU PUT PATIENTS' HEALTH FIRST,  
REMARKABLE THINGS FOLLOW



State-of-the-art manufacturing facilities

OUR RANGE OF PRODUCTS INCLUDE

SINGLE USE CIRCULAR STAPLER FOR RECTAL PROLAPSE AND HAEMORRHOIDS	SINGLE USE RELOADABLE LINEAR STAPLER	SINGLE USE CURVED INTRALUMINAL CIRCULAR STAPLER	CPH34HV	CS COMPACT™	SINGLE USE TRANSVERSE CUTTING LINEAR STAPLER
--	--------------------------------------	---	---------	-------------	--



FOR MORE INFORMATION OR TO VIEW OUR ENTIRE PRODUCT RANGE PLEASE VISIT OUR WEBSITE AT [WWW.FRANKENMAN.COM](http://WWW.FRANKENMAN.COM)

CE 0123



FRANKENMAN INTERNATIONAL LTD.  
SUITE B, 13/F, WING TAT COMMERCIAL BUILDING  
121-125 WING LOK STREET, SHEUNG WAN, HONG KONG  
TELEPHONE: +852 3106 3035 | FAX: +852 3585 0519  
E-MAIL: [info@frankenman.com](mailto:info@frankenman.com)

[www.frankenman.com](http://www.frankenman.com)



**FRANKENMAN**

INFORMED DESIGN.  
SCRUPULOUSLY APPLIED.  
METICULOUSLY PRODUCED.



SINGLE USE RELOADABLE  
LINEAR CUTTER STAPLER  
AND RELOADS

# A RANGE OF PRODUCTS TO MEET THE NEEDS OF MODERN SURGICAL PRACTICE

The range of Linear Cutter Staplers focuses on reliability and functionality. Designed to provide accurate and reliable staple formation and tissue resection... time and time again.

At Frankenman, we believe we're not so much in the business of making instruments, but in the business of making surgeries more successful. That's why we listen to doctors who've put our instruments to use. We read their reports. We invite them to think with us, to brainstorm, to help us develop new and innovative ideas. That way we know we will be able to deliver the highest quality products to the medical profession both today and way into the future.

## THE LC RANGE - KEY FEATURES

1

SIDE VIEW - LC60/80



1

TISSUE RETAINING PIN

2

CLEAR GRADUATIONS AND INDICATION OF 'CUTTING ZONE' ON BOTH SIDES OF STAPLE HOUSING

3

ALL RE-LOADABLE CARTRIDGES HAVE SAFETY LOCK TO PREVENT FIRING OVER EMPTY CARTRIDGE

6

MODERN ERGONOMIC DESIGN DEVELOPED WITH INTERNATIONALLY RENOWNED BRITISH INDUSTRIAL DESIGNERS

7

UNIQUE FIRING 'SAFETY LOCK' TO PREVENT INADVERTENT PRE-FIRING DURING PREPARATION AND POSITIONING

4

UNIQUE 90 DEGREE RELOAD 'SHOULDER' TO ACT AS A POSITIVE TISSUE LIMITER TO AVOID CONCERTINA EFFECT OVER RE-LOAD AND RISK OF PRE-CUT INTERNALLY

SIDE VIEW-LC100



5

CARTRIDGE REMOVAL MADE EASIER. INDENTATION TO EXPOSE CARTRIDGE "LEGS" FACILITATES EASY REMOVAL WHEN CARTRIDGE IS WEDGED BY TISSUE OR FLUID BUILD UP



## PRODUCT SPECIFICATIONS

PRODUCT CODE				CLOSED STAPLE SIZE	STAPLING LINE LENGTH	STAPLE QUANTITY
LINEAR CUTTER STAPLER		CARTRIDGE RELOADS				
LC60G	GREEN	LCR60G	GREEN	2.0MM	61MM	60
LC60B	BLUE	LCR60B	BLUE	1.5MM	61MM	60
LC60W	WHITE	LCR60W	WHITE	1.0MM	61MM	60
LC80G	GREEN	LCR80G	GREEN	2.0MM	81MM	80
LC80B	BLUE	LCR80B	BLUE	1.5MM	81MM	80
LC100G	GREEN	LCR100G	GREEN	2.0MM	101MM	100
LC100B	BLUE	LCR100B	BLUE	1.5MM	101MM	100

## KEY FEATURES

- CAM GEAR DESIGN: ENSURES CONSISTENT STAPLE FORMATION HELPING PREVENT BLEEDING.
- EMPTY CARTRIDGE SAFETY LOCK: CUTTING BLADE WILL ONLY ACTIVATE WITH FULL CARTRIDGE.
- RELOADABLE COLOUR CODED CARTRIDGES.